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MENU OF AVAILABLE HEART SOUNDS

<u>LOCATION</u>	<u>HEART SOUND</u>	<u>COMMENT</u>
Base Right	Venous Hum	This continuous murmur may be found in children aged 3 to 6 years. It occurs as a result of turbulence in the jugular venous system and is only heard when the child is in the upright position. The rate is about 96 beats per minute.
Base Right	Aortic Stenosis	This systolic murmur is loudest over the ascending aorta. Duration and intensity vary with the severity of stenosis. An ejection click may be heard..
Base Left	Split S2	Sounds are rather normal. Degree of splitting increases with inspiration and decreases with expiration. Wide split suggests prolonged RV ejection or shortened LV ejection; narrow split suggests early closure of pulmonary valve.
Base Left	Pulmonary Stenosis	This systolic murmur is normally loudest over the main pulmonary artery. Duration and intensity vary with the severity of stenosis.
Base Left	Systolic Fixed S2	Fixed S2 does not change width during respiration. The absence of split S2 usually indicates a condition that prolongs RV ejection time or shortens LV ejection. Conditions include volume/pressure overload and RBBB.
Apex	One Year Heart	Normal heart sounds heard in a 12 month old infant. The rate is 120 beats per minute.
Apex	Six Year Heart	Normal heart sounds heard in a 6 year old child. The rate is 84 beats per minute.

MENU OF AVAILABLE HEART SOUNDS (CONTINUED)

<u>LOCATION</u>	<u>HEART SOUND</u>	<u>COMMENT</u>
Apex	Stills Murmur	The vibratory murmur may be found in children between ages 3 to 6 years and sounds like "twanging string", or squeaking/buzzing at a low frequency. During inspiration, murmurs increase on the right side and decrease on the left.
Apex	Mitral Valve Regurgitation	This systolic murmur is produced by lesions, more often rheumatic than congenital in origin. Chest films may show pulmonary vein congestion, pulmonary edema, or an enlargement of the left atrium/ventricle.
Apex	Split S1	S1 relates to closure of mitral and tricuspid valves. Not common in normal children, and may indicate RBBB or other anomalies. Sound can be confused with an ejection click on S4.

MENU OF AVAILABLE LUNG SOUNDS

<u>LOCATION</u>	<u>LUNG SOUNDS</u>	<u>COMMENT</u>
Trachea	Normal Infant	Expiration sounds are louder, have a higher pitch, and are of longer duration than during inspiration. The silent period or pause following expiration is longer than the one between expiration and inspiration.
Trachea	Normal Child	Expiration sounds are louder. Have a higher pitch, and are of longer duration than during inspiration. The silent period or pause following expiration is longer than the one between expiration and inspiration.
Trachea	Stridor Sounds	Patient has marked respiratory distress, and a narrow aperture between the vocal cords that produces a high pitched tone during both inspiration and expiration. During the end of expiration, there is an abrupt drop in pitch.
Upper Anterior (2 sites) and Lower Anterior (2 sites)	Wheezing Sounds	These musical wheezing sounds are often heard in asthma patients. During inspiration, the wheeze is slightly higher in pitch than during expiration. Wheezing in asthmatics is often present in either one or both phases of respiration.
Upper Posterior (2 sites) and Lower Posterior (2 sites)	Crackles, Ronchi	Coarse crackles are present during both inspiration and expiration. There are also some very low pitched repetitive sounds that are ronchi. High pitched squeaks are also audible against a background of bronchial breath sounds.

DESCRIPTION OF AVAILABLE SOUNDS LOCATED
ON THE PEDIATRIC VIRTUAL STETHOSCOPE

BASE RIGHT	TRACHEA
AORTIC STENOSIS	NORMAL INFANT
AORTIC STENOSIS	NORMAL CHILD
VENOUS HUM	STRIDOR SOUNDS
VENOUS HUM	STRIDOR SOUNDS
VENOUS HUM	STRIDOR SOUNDS
BASE LEFT	UPPER ANTERIOR
SPLIT S2	WHEEZING SOUNDS
SYSTOLIC FIXED S2	WHEEZING SOUNDS
SYSTOLIC FIXED S2	WHEEZING SOUNDS
PULMONIC STENOSIS	WHEEZING SOUNDS
PULMONIC STENOSIS	WHEEZING SOUNDS
LLSB	LOWER ANTERIOR
PULMONARY STENOSIS	WHEEZING SOUNDS
PULMONARY STENOSIS	WHEEZING SOUNDS
SPLIT S1	WHEEZING SOUNDS
SPLIT S1	WHEEZING SOUNDS
SPLIT S1	WHEEZING SOUNDS
APEX	POSTERIOR
1 YEAR HEART	RONCHI CRACKLES
6 YEAR HEART	RONCHI CRACKLES
STILLS MURMUR	RONCHI CRACKLES
SPLIT S1	RONCHI CRACKLES
MITRAL VALVE REGURG.	RONCHI CRACKLES